

FOR PARTNER CLINICS USING 3SHAPE TRIOS

Send TRIOS scans to ASO Hawaii *via EasyRx.*

Export STL from TRIOS, upload to EasyRx, done -
secure, HIPAA-compliant, with full case history.

1 One-time setup

Free EasyRx account
invited by ASO

2 Export & upload

TRIOS STL export,
then drag into EasyRx

3 Full case history

HIPAA-compliant,
no file size limits

Why EasyRx instead of email?

- **HIPAA-compliant** - patient info stays encrypted end-to-end
- **No file size limits** - TRIOS STL files upload without truncation
- **Case tracking** - see status from submission to delivery
- **Prescription together** - Rx form and STL files stay linked

EMAIL STILL WORKS AS A BACKUP

If EasyRx is temporarily unavailable, email STL files to aso-digital@outlook.com. See Page 3.

ONE-TIME SETUP · ~5 MINUTES

Get EasyRx + export from TRIOS

Two short setups - once done, every case takes about 2 minutes.

1

Request EasyRx access from ASO

Visit asohawaii.com/contact and fill out Request an EasyRx Connection. ASO sends you an EasyRx invitation email within 1 business day - it's free for you.

2

Accept the invitation and set your password

Click the link in the EasyRx email, set your password, and you're ready to log in at app.easyrxortho.com. No license fees for clinic accounts.

3

In TRIOS Dental Manager: enable STL export

Open Dental Manager. Make sure the "TRIOS Scan STL Export" option is enabled on your dongle. Ask your 3Shape reseller if it isn't (one-time dongle setting).

4

Export the scan as STL

In Dental Manager, right-click the TRIOS order, choose Advanced → Export Scan. Pick STL format, choose Z-axis, select a folder, and save.

HAVE 3SHAPE ORTHOANALYZER?

You can set up auto-export to a folder so STLs save automatically. Contact ASO for setup help.

NAMING TIP

Rename STL files to include patient name (e.g. "Smith_John_Upper.stl") before uploading.

EVERY CASE · ~2 MINUTES

Submit your TRIOS case

After exporting STL from TRIOS, follow these 3 steps in EasyRx.

5

Log in to EasyRx

Go to app.easyrxortho.com and sign in. Click Create Rx.

6

Enter patient + Rx

Add patient name, select appliance type, due date, and any notes.

7

Upload STL → Send

Drag your STL files into the Rx, then click Submit to send to ASO.

Email fallback (only if EasyRx is down)

If you cannot access EasyRx temporarily, email STL files with these details:

TO: aso-digital@outlook.com

SUBJECT: TRIOS Case - [Patient Name] - [Appliance Type]

INCLUDE IN BODY:

- Patient name, DOB, doctor name
- Appliance type, due date, special instructions
- Attach both upper/lower STL files (rename with patient name)

ASO HAWAII · CASE QUESTIONS

Tel · 808-957-0111

Email · aso-digital@outlook.com

Hours · Mon-Fri, 8:00 AM - 4:30 PM HST

EASYPX & 3SHAPE SUPPORT

EasyRx · support@easyrxortho.com

EasyRx phone · 1-888-340-3751

3Shape · support.3shape.com