

FOR PARTNER CLINICS USING DEXIS IS 3800 / 3800W

Send scans to ASO Hawaii *via IS Connect.*

Using your DEXIS IS 3800 / 3800W intraoral scanner -
one connection, one-click transfer, straight into our queue.

1

One-time setup

Request connection
in IS Connect portal

2

Direct from scan

IS ScanFlow sends
cases with one click

3

Auto-received

ASO polls IS Connect
every 5 minutes

Before you start

- **A DEXIS SSO account** - the same account you use for IS ScanFlow
- **IS Connect access** - free - included with your DEXIS IS 3800 / 3800W
- **About 2 minutes** - for the one-time partnership connection

IS CONNECT CLOUD OR CLASSIC?

DEXIS launched IS Connect Cloud in 2026 - both versions work with ASO.
The steps below apply to either version. UI labels may vary slightly.

ONE-TIME SETUP · ~2 MINUTES

Connect your clinic to ASO International Hawaii

You only need to do this once per clinic account.

1

Sign in to DEXIS IS Connect

Open IS Connect at dentalconnect.dexis.com or launch it from the IS ScanFlow toolbar. Sign in with your DEXIS SSO account.

2

Open Partners

Click the Partners icon in the top toolbar (the icon with two people).

3

Request connection to our lab

Click Add Partner / Invite Lab and enter our email address below. DEXIS will send us a connection request.

aso-digital@outlook.com

4

Wait for approval

We accept within 1 business day. You'll see ASO International Hawaii in your Partners list once the connection is active. You're then ready to send cases.

WHAT HAPPENS NEXT

1. ASO receives your connection request by email (usually accepted same day)
2. You'll see "ASO International Hawaii" added to your Partners list in IS Connect
3. You're ready to start sending cases - see Page 3

EVERY CASE · ~30 SECONDS

Send a case to ASO

After scanning with IS ScanFlow, follow these 4 steps.

1

Finish scanning

Complete your scan in IS ScanFlow. Confirm quality on the preview screen and click Finalize or Save.

2

Choose Send to Lab

From IS ScanFlow, click the Send icon or choose Transfer to IS Connect from the workflow menu.

3

Select ASO International Hawaii

Pick ASO International Hawaii from the Partners dropdown. Your saved prescription form will auto-load.

4

Fill Rx details -> Send

Enter appliance type, due date, and notes. Click Send. Status changes to Transferred - your case is now in our queue.

Case status in IS Connect

Transferred

->

Consulted

->

In Production

->

Shipped

Track any case in IS Connect - your case moves through these stages until delivery.

ASO HAWAII · CASE QUESTIONS

Tel · 808-957-0111

Email · aso-digital@outlook.com

Hours · Mon-Fri, 8:00 AM - 4:30 PM HST

DEXIS · TECHNICAL SUPPORT

Web · dexas.com/support

Phone · 1-888-883-3947

Portal · dentalconnect.dexas.com