

FOR PARTNER CLINICS USING MEDIT

Send cases to ASO Hawaii *in 4 taps.*

Using Medit Link with your i700 / i900 intraoral scanner —
no email, no export, no paperwork.

1 One-time setup

Partnership connection
in under a minute

2 Direct transfer

From Medit Link
straight to ASO

3 Auto-received

Polled every 5 minutes
to our EasyRx queue

Before you start

- **A Medit Link account** — free — sign up at meditlink.com if you haven't
- **Your Medit intraoral scanner** — connected and working (i700 / i900 recommended)
- **About 2 minutes** — for the one-time partnership setup

★ ALREADY A PARTNER?

Skip to Page 3 — the 4-step send workflow is all you need each time.

ONE-TIME SETUP · ~2 MINUTES

Connect your clinic to ASO International Hawaii, Inc.

You only need to do this once per clinic account.

1

Open Medit Link

Launch the Medit Link app (desktop or web at meditlink.com) and sign in with your clinic account.

2

Go to Partners → Find Partners

In the left sidebar, click the Partners icon, then select Find Partners at the top.

3

Search for our lab

Type the following into the search bar — exactly as shown:

ASO International Hawaii

4

Click Request Partnership

Find our lab in the results (look for the ASO logo — a gold tooth on navy background) and click Request Partnership. We'll accept within 1 business day.

WHAT HAPPENS NEXT

1. ASO Hawaii accepts your partnership request (usually same day, within 1 business day)
2. You'll see a notification in Medit Link and ASO will appear in your Partners list
3. You're ready to start sending cases — see Page 3

EVERY CASE · ~30 SECONDS

Send a case to ASO

Finish your scan, then follow these 4 steps — nothing else required.

1

Open the scanned case

In Medit Link, go to Case Box and select the case you want to send. Confirm scan quality on the final review screen.

2

Click Order → Send to Lab

From the case detail page, click the Order button at the top right, then choose Send to Lab.

3

Select ASO International Hawaii

In the lab dropdown, pick "ASO International Hawaii, Inc." Your EasyRx prescription form will auto-populate.

4

Fill Rx details → Send

Enter appliance type, due date, and notes. Click Send. Your case is now in our queue and you'll see "Waiting to be accepted."

After submission



When your package arrives, click Confirm Delivery in Medit Link to close the case.

ASO HAWAII · CASE QUESTIONS

Tel · 808-957-0111

Email · aso-digital@outlook.com

Hours · Mon–Fri, 8:00 AM – 4:30 PM HST

MEDIT LINK · TECHNICAL SUPPORT

Web · support.medit.com

Email · support@medit.com

In-app · Help → Contact Us