

FOR PARTNER CLINICS USING PRIMESCAN

# Send Primescan scans to ASO *via DS Core.*

Primescan uploads to DS Core automatically -  
share with ASO and we receive your case instantly.

**1**

## Partnership setup

One-time connection  
request in DS Core

**2**

## Scan and share

Primescan -> DS Core  
-> ASO in 4 clicks

**3**

## Auto-received

ASO gets it instantly  
via DS Core

## Before you start

- **A DS Core account** - the same account you use to sync your Primescan
- **Primescan connected to DS Core** - check Settings -> Connect in Primescan Software
- **About 2 minutes** - for the one-time partnership with ASO

### WHAT IS DS CORE?

DS Core is Dentsply Sirona's cloud platform. Your Primescan uploads scans there automatically.  
Access it at [dscore.dentsplysirona.com](https://dscore.dentsplysirona.com) or from the DS Core icon in Primescan Software.

ONE-TIME SETUP · ~2 MINUTES

# Connect your clinic to ASO International Hawaii

You only need to do this once per clinic account.

1

## Sign in to DS Core

Open [dscore.dentsplysirona.com](https://dscore.dentsplysirona.com) or launch DS Core from Primescan Software. Sign in with your DS Core account.

2

## Go to Share -> Partners

In the left sidebar, click Share (the share icon). Choose the Partners tab to manage your lab connections.

3

## Search for our lab

Click Add Partner and type the following into the search bar:

**ASO International Hawaii**

4

## Send connection request

Select our lab from the results and click Request Partnership. We'll accept within 1 business day - you'll get an email when we do.

### WHAT HAPPENS NEXT

1. ASO Hawaii accepts your connection request (usually same day)
2. "ASO International Hawaii, Inc." appears in your Partners list in DS Core
3. You're ready to start sending cases - see Page 3

EVERY CASE · ~30 SECONDS

# Send a case to ASO

After your Primescan uploads to DS Core, share in 4 clicks.

1

## Open the patient in DS Core

In DS Core, go to Patients and select the patient. Your Primescan scan is already uploaded.

2

## Click Share -> Send to Lab

On the patient screen, click the Share icon. Choose Send to Lab from the dropdown.

3

## Select ASO International Hawaii

Pick ASO International Hawaii from your Partners list. Your prescription form will load automatically.

4

## Fill Rx and send

Enter appliance type, due date, and notes. Click Send. Status will show Sent - you are done.

## Case status in DS Core

Sent

->

Accepted

->

In Production

->

Shipped

Track any case in DS Core - your case moves through these stages until delivery.

### ASO HAWAII · CASE QUESTIONS

Tel · 808-957-0111

Email · aso-digital@outlook.com

Hours · Mon-Fri, 8:00 AM - 4:30 PM HST

### DS CORE · TECHNICAL SUPPORT

Web · support.dentsplysirona.com

Portal · dscore.dentsplysirona.com

Phone · 1-800-659-5977