

FOR PARTNER CLINICS USING SHINING 3D SCANNERS

Send Shining 3D scans to ASO *via Dental Cloud.*

Aoralscan and DS-EX Pro scanners upload to Dental Cloud - share with ASO as a partner for direct case delivery.

1

Partnership setup

One-time connection
request in Dental Cloud

2

Scan and share

Aoralscan -> Cloud
-> ASO in 4 clicks

3

PLY color retained

Vertex color data
preserved end-to-end

Before you start

- **A Dental Cloud account** - free - comes with your Shining 3D scanner
- **Your scanner linked to the cloud** - Aoralscan 3 / Elite, DS-EX Pro, or newer
- **About 2 minutes** - for the one-time partnership with ASO

WHAT IS SHINING 3D DENTAL CLOUD?

Dental Cloud is Shining 3D's case-management and transfer platform for intraoral scanners.
Access it at dentalcloud.shining3d.com or through Aoralscan / DS-EX Pro software.

ONE-TIME SETUP · ~2 MINUTES

Connect your clinic to ASO INTERNATIONAL HAWAII inc

You only need to do this once per clinic account.

1

Sign in to Dental Cloud

Open dentalcloud.shining3d.com in your browser, or launch it from your Aoralscan / DS-EX Pro software. Sign in with your clinic account.

2

Go to Partners

In the left sidebar, click the Partner icon. You'll see your current lab connections (empty on first visit).

3

Search for our lab

Click Add Partner or Find Partner. In the search field, type:

ASO INTERNATIONAL HAWAII inc

4

Send connection request

Select our lab from the results (Hawaii, United States) and click Apply / Connect. We'll approve within 1 business day.

WHAT HAPPENS NEXT

1. ASO Hawaii receives your connection request (usually accepted same day)
2. "ASO INTERNATIONAL HAWAII inc" appears in your Partners list in Dental Cloud
3. You're ready to start sending cases - see Page 3

EVERY CASE · ~30 SECONDS

Send a case to ASO

After your Aoralscan uploads to Dental Cloud, share in 4 clicks.

1

Open the order in Dental Cloud

Go to Orders in Dental Cloud and select the case. Your Aoralscan upload is already attached.

2

Click Send / Share

On the order screen, click Send to Lab or the share icon to start the transfer workflow.

3

Select ASO INTERNATIONAL HAWAII

Pick ASO INTERNATIONAL HAWAII inc from your Partners list. Fill in the prescription details.

4

Confirm and send

Review appliance type, due date, and any notes. Click Send. Status updates to Submitted - you are done.

Case status in Dental Cloud

Submitted

->

Accepted

->

In Production

->

Shipped

Track any case in Dental Cloud - your case moves through these stages until delivery.

ASO HAWAII · CASE QUESTIONS

Tel · 808-957-0111

Email · aso-digital@outlook.com

Hours · Mon-Fri, 8:00 AM - 4:30 PM HST

SHINING 3D · TECHNICAL SUPPORT

Web · dental.shining3d.com

Cloud · dentalcloud.shining3d.com

Email · dental_support@shining3d.com