

iTero + EasyRx Setup Guide

For practices using the iTero intraoral scanner to submit cases to ASO Hawaii

What this guide covers

Once set up, your iTero scans will be **automatically attached** to every EasyRx prescription you submit to ASO Hawaii. No more downloading STL files, no more email attachments, no more lost scans.

Setup takes about **5 minutes**, and only needs to be done once per practice.

Before You Begin

You will need the following:

- An **EasyRx account** connected to ASO International Hawaii (if you don't have one yet, contact ASO Hawaii at aso-digital@outlook.com to receive a free invitation)
- Your **MyiTero.com login credentials** (email address and password)
- **Admin access** to your EasyRx account (the integration can only be set up by an admin user)
- A web browser: **Chrome or Firefox** is recommended by EasyRx

Cost to Your Practice

Free. ASO International Hawaii covers 100% of the EasyRx subscription cost for your practice. You will not be asked for a credit card or any payment information at any point.

Part 1 — Enable the iTero Integration in EasyRx

Complete these 5 steps once, as the Admin user of your EasyRx account.

1 Log in to EasyRx as an Admin

Go to **easyrxortho.com** and sign in with your Admin credentials. If you are not sure whether your account has Admin rights, ask your practice manager or contact EasyRx support at 1-888-340-3751.

2 Open Account Settings

In the **top-right corner**, click on your practice name (shown in blue). In the dropdown menu, click **Account Settings**.

3 Go to 3rd Party Integrations

On the left side of the screen, click **3rd Party Integrations**.

If you do not see this option, contact EasyRx support — they can confirm your admin email and enable the option for you.

4 Add iTero Credentials

In the **iTero Credentials** section, click **Add**. You will be redirected to an authentication page where you will enter your **MyiTero.com email and password**.

Follow the prompts to complete authentication and return to EasyRx.

5 Pair Your Office

After authentication, you will see a list of iTero offices associated with your MyiTero account. Click **Pair Office** next to each office you want to link to this EasyRx account.

Once paired, all scans from those offices will appear automatically in your iTero Recent Scans dashboard in EasyRx.

v1 vs v2 — What you need to know

EasyRx recently upgraded from v1 to **v2** of the iTero integration for all U.S., Canada, and Mexico users. v2 authenticates directly with MyiTero.com (not MyAlignTech). If you completed setup a long time ago and still see a "v1" label next to your scans, follow the steps above to re-authenticate with v2.

Part 2 — Submit a Case with an iTero Scan

Once the integration is set up, submitting a case takes just a few clicks. There are two equally easy methods — choose whichever is more convenient.

Method A — Start from the iTero Dashboard (Recommended)

This is the fastest way because the scan is linked automatically.

1 Take the scan on your iTero scanner

Use either **iRecord** or **iRecord + Invisalign** scan type. These scan types upload to your MyiTero account within about 30 minutes and work seamlessly with the EasyRx integration.

Note: You do NOT need to select a lab on the iTero scanner itself. The EasyRx integration handles routing to ASO Hawaii.

2 Open the iTero — Recent Scans dashboard in EasyRx

From the EasyRx main dashboard, click **iTero — Recent Scans**. You will see a list of all scans from your paired offices, with the most recent scans at the top.

If your scan is not yet showing, click **Check for New Scans** at the top.

3 Click Create Rx

Find your patient's scan in the list and click **Create Rx** on the right-hand side.

- If the patient already exists in EasyRx, the Create Rx screen opens directly.
- If the patient is new, EasyRx will prompt you to create the patient record first.

4 Fill out the prescription and submit

Complete the appliance selection, design instructions, and due date. The iTero Case Order ID is already filled in, so your scan will be attached automatically when you click **Submit**.

Scan types and turnaround

iRecord / iRecord + Invisalign: scans typically available within 30 minutes, auto-attach is immediate.

iCast scans: require 2–5 business days of processing by the MyAlignTech team before EasyRx can attach them. For iCast, please allow at least 10 calendar days before the due date to be safe.

Method B — Start from a New Rx, Then Attach the Scan

If you prefer to start from the prescription form, you can attach the iTero scan from the Uploads tab.

1 Create a new Rx as usual

From the EasyRx dashboard, click **New Rx** and select the patient. Fill out the appliance and design sections as you normally would.

2 Go to the Uploads tab

Inside the Rx, click the **Uploads** tab. You will see an **iTero Case Order ID** field along with a **Lookup and Insert iTero Case** button.

3 Look up the scan

Click **Lookup and Insert iTero Case**. A window will open showing your recent iTero scans. Find the correct patient and click **Attach Scan**.

*Alternatively, if you already know the Case Order ID, type it into the iTero Case Order ID field and click **Validate**.*

4 Submit the prescription

Once validated, the scan will auto-attach as an STL file when you submit. The attachment process usually takes a few minutes.

Part 3 — Verifying Your Submission

After you submit, you can confirm that your scan was attached successfully.

- Prescriptions with iTero scans are displayed in **green** in the Rx list.
- To view attached files, open the Rx and go to the **Uploads** tab.
- If the auto-attach is still in progress, EasyRx will display a message. Refresh the page after a few minutes to see the STL files.
- ASO Hawaii will begin processing your case as soon as the scan is received — no additional email or file transfer is needed.

HIPAA Compliance

All files uploaded through EasyRx servers are backed up in accordance with HIPAA and ADA requirements.

Troubleshooting

Issue	What to do
"Unable to validate this Case Order ID"	The Case Order ID was entered incorrectly or has extra spaces. Double-check the ID in your MyiTero account and re-enter it. Case IDs are always numeric — there are no letters.
Scan is not showing up in EasyRx	• Confirm the scan is marked "Completed" in MyiTero.com. • Click "Check for New Scans" in the iTero Recent Scans dashboard. • Verify that the office was successfully paired (Part 1, Step 5). • For iCast scans, wait 2–5 business days for MyAlignTech processing.
Scans are labeled "v1" in the dashboard	These scans were attached using the older v1 integration. They will continue to work for a limited time, but you should re-authenticate using the v2 steps in Part 1 of this guide.
Can't see "3rd Party Integrations" in Account Settings	You are likely not logged in as an Admin. Contact your EasyRx account admin, or call EasyRx support at 1-888-340-3751 to confirm your role.
Need to change the base type of a scan	For iRecord / iRecord + Invisalign scans, ASO Hawaii can re-export with OrthoCAD to add a low-profile base. For iCast scans, the base is already applied and cannot be changed after the fact.

Support Contacts

Contact	Details
ASO International Hawaii	Phone: 808-957-0111 Email: aso-digital@outlook.com Address: 1441 Kapiolani Blvd #1112, Honolulu, HI 96814 For: case questions, appliance design, EasyRx invitation requests
EasyRx Support	Phone: 1-888-340-3751 Email: support@easyrxcloud.com For: login issues, integration setup, account access
iTero / MyiTero Support	Contact your Align Technology representative or visit myitero.com For: scanner issues, MyiTero login recovery, scan processing